MORVABON

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manufacturer

Pardis Morvabon Pharmed Co.:

Head Office: Unit 154, 7th floor, c block,
Borjsaz Bd, Azadi Street, Tehran, Province, Iran

Tel: 0098 21 6658 1838

Fax: 0098 21 6690 8819

Factory: Unit 135,136,139, 2nd floor, Seraj Bd.,
20thkm of Damavand Rd., Tehran, Iran

Tel: 0098 21 7625 1702

Fax: 0098 21 7648 1457

MEDICAL
EQUIPMENT
DEPARTMENT

Glass Ionomer Luting Cement

Radiopaque self cure glass ionomer luting and lining cement

Morvabon Light Cure Glass Ionomer is a radio opaque cement that contains fluoride and has high compressive strength and adhesive strength to dentin. It is a suitable product for use as a base and liner under composite, amalgam or porcelain.

Composition:

Fluoroaluminosilicate , Polymer , initiator

Application:

Can be used as a base and liner
Direct repair of small damages

Advantages:

- ▶ High radiopacity
- ▶ Easy handling
- ▶ Match the color of tooth dentin
- ▶ High compressive strength
- ▶ Contains fluoride

Precautions:

- ▶ Carefully read and understand instruction before using.
- ▶ Avoid contact with oral soft tissue, eyes or skin.
- ▶ in case of consumption drink copious amount of water.do not induce vomiting. seek medical advice.
- ▶ if accidental contact with soft tissues, eyes or skin occurs, immediately flush with plenty of water.
- ▶ Use eye protection, mask and gloves when using.
- ▶ It is necessary to use a rubber band while working with this product.
- ▶ Do not use for patients who are highly allergic to the components of this product.
- ▶ For dental use only.

Warning:

- ▶ Store at 5-25 °C.
 - ▶ Keep away from moisture and direct sunlight.
 - ▶ Close the lid of the product after use to prevent absorption.
- Do not use after expiration date.
Keep out of reach of children.

Contraindication:

- ▶ This product can be used by all people who are not allergic to its active ingredient.
- ▶ Do not use in pulp capping technique.

Instructions for use:

1. Prepare the cavity in the usual way.
2. Rinse and dry the cavity.
3. Dispense the material on the pad.

Note: Avoid direct light exposure to the dispensed material.

4. Apply the product by the applicator to obtain a first coating of no more than 0.3 to 0.5 mm thickness.
5. Remove the material that inadvertently applied over enamel margins.
6. Cure for 20 seconds (30 seconds with devices with lower light intensity).

Note: This product is intended to be cured by exposure to a halogen or LED light with a minimum intensity of 400 mW/cm² in the 400-500 nm range.

7. If intended, apply the second layer and cure for 30 seconds. The second layer can be up to 1.5 mm thick.

 Consult instructions for use	 Indicates the batch code	 Indicates the date after which the medical device not to be used.	 Indicates a medical device that needs protection from light sources.
 Do not use if Package is Damaged	 Indicates the medical device manufacturer (Name & Address)	 Indicates a medical device that has not been subjected to a sterilization process.	 Indicates the device intended for use by professional person only.
			 Indicates a medical device that needs to be protected from moisture.

Produced by Morvabon Co.
www.morvabon.com
info@morvabon.com
Made in IRAN